Get the Blood Drop - Palm

ALTERNATIVE SITE TESTING (AST) PALM Use the clear AST endcap. The AST endcap is not recommended for pediatric population. (Puncture depth = 3.33 mm, non-adjustable)

- **1.** Insert the lancet as described in *Prepare* the Lancing Device.
- 2. Attach and lock the clear AST endcap. Twist the white band until you hear or feel a click, and the white and black raised marks line up as shown.



CAUTION: Do not use the clear AST endcap on your fingertips.

3. Select a puncture site from a fleshy area on the palm. Avoid veins, moles, bones, and tendons.

- 4. Pull back the re-set handle to prepare the lancing device. The release button is raised when the lancing device is ready.
- 5. Press the lancing device firmly against the puncture site, then press the release button UNTIL THE DEVICE CLICKS.




- (If the device does not puncture, reset by repeating Step 4.)
- **6.** Do not take the device away from your palm. Maintain steady pressure until a small, round blood drop forms under the clear endcap.



7. Lift the lancing device straight up and away from the skin without smearing the blood drop.

Test immediately.



CAUTION: Do not test the blood sample from palm if you get: Smeared blood

- Runny blood
- Clotted blood

Clear fluid mixed with the blood

9. Wipe the puncture site with a clean, dry tissue and maintain pressure on the site until the bleeding stops.

See: Eject / Dispose of the Lancet

Cleaning and Disinfection

WARNING: Potential Biohazard

If your lancing device is being operated by a second person who is providing testing assistance to you, the second person should wear gloves and the lancing device should be disinfected prior to use to reduce the risk of transmitting blood-borne pathogens.

Clean and disinfect your lancing device after each use. Use only Super Sani-Cloth germicidal wipes containing 55% isopropyl alcohol and 0.5% ammonium chloride. The wipes are proven safe to use with the lancing device.

Cleaning involves the removal of visible dirt and debris, but does not reduce the risk for transmission of infectious diseases.

Disinfecting (if performed properly) reduces the risk of transmitting infectious diseases.

Super Sani-Cloth germicidal wipes are available for purchase online. For more information, call Customer Service at 1-800-348-8100, available Monday through Sunday, from 8:00 am through 12:00 Midnight, Eastern Time.

CAUTION: Do not use bleach (sodium hypochlorite).

For more information:


Centers for Disease Control and Prevention. CDC clinical reminder: use of fingerstick devices on more than one person poses risk for transmitting bloodborne pathogens. US Department of Health and Human Services; 8/23/2010. http://www.cdc.gov/ injectionsafety/Fingerstick-DevicesBGM.html US Food and Drug Administration. Use of fingerstick devices on more than one person poses risk for transmitting bloodborne pathogens: initial communication. US Department of Health and Human Services; update 11/29/2010. http://wayback. archive-it.org/7993/20170111013014/http://www. fda.gov/MedicalDevices/Safety/AlertsandNotices/ ucm224025.htm

Signs of Deterioration: Your lancing device has been tested for 3,650 cycles of the cleaning and disinfection procedure (5 cycles per day for 2 years). The lancing device has been demonstrated to withstand 2 years of cleaning and disinfection without damage. However, if your device malfunctions for any reason, contact Customer Service at 1-800-348-8100 available Monday through Sunday, from 8:00 am through 12:00 Midnight, Eastern Time to troubleshoot and provide a solution.

Cleaning Your Lancing Device

Supplies Needed for Cleaning:

- Super Sani-Cloth germicidal wipes
- . Remove the endcap and eject the lancet. if presen[.]
- **2.** Clean the lancing device body with Super Sani-Cloth germicidal wipes until all soil is removed.



- Clean the endcap inside and outside until all soil is removed.
- Visually inspect the lancing device body and endcap for any residual soil. Repeat steps 2 and 3 if any soil is still visible.

See image on next page



5. Allow to air dry.

Disinfecting Your Lancing Device

Supplies Needed for Disinfectina: Super Sani-Cloth germicidal wipes

- Timing device
- Before disinfecting, clean your lancing device as described in Cleaning Your Lancing Device.
- 2. Get a new Super Sani-Cloth germicidal
- 3. Set timer for 2 minutes and start timing.
- LANCING DEVICE BODY: Wipe the body continuously. You must wipe all surfaces and keep them wet for 2 minutes.

- 5. Set timer for 2 minutes and start timing
- 6. ENDCAP: Wipe the endcap continuously. You must wipe the inside and outside and keep them wet for 2 minutes.
- Allow lancing device body and endcap to air dry.

Storage

Store your lancing device in your meter carrying case.

Customer Service

Warranty: Your lancing device is warranted for 2 years.

If you need assistance or have questions about blood sugar testing, please contact Customer Service at 1-800-348-8100, available Monday through Sunday, from 8:00 am through 12:00 Midnight, Eastern Time.

MICROLET NEXT lancing device made in Poland.

UDI = Unique device identifier

Ascensia Diabetes Care US Inc. 5 Wood Hollow Road Parsippany, NJ 07054 Tel: 1-800-348-8100 8:00 AM to 12:00 AM (ET)

Español

Para instrucciones en Español comuníquese a la línea del servicio al cliente 1-800-348-8100 o visite la página www.diabetes.ascensia.com

ASCENSIA Diabetes Care

For information on patents and related licenses, see: www.patents.ascensia.com

Ascensia, the Ascensia Diabetes Care logo, and Microlet are trademarks and/or registered trademarks of Ascensia Diabetes Care Holdings AG, All other trademarks are properties of their respective owners and are used solely for informative purposes. No relationship or endorsement should be inferred or


© 2022 Ascensia Diabetes Care Holdings AG. All rights reserved.





INSTRUCTIONAL INSERT

- Lancet ejector for safe lancet disposal
- Black/white adjustable endcap for fingertip testing only.
- Clear AST endcap for alternative site testing only.
- The MICROLET®NEXT lancing device and MICROLET[®] lancets are used to draw capillary blood from the fingertip or palm for blood glucose testing.





区

WARNING

- Ask your health care professional if Alternative Site Testing (AST) is right for you.
- Do not calibrate a continuous glucose monitoring device from an AST result.
- Do not calculate an insulin dose based on an AST result.

Do not use AST under the following conditions:


- If you think your blood sugar is low (hypoglycemia).
- When blood sugar is changing rapidly (after a meal, insulin dose, or exercise).
- If you are unable to feel symptoms of low blood sugar (hypoglycemic unawareness).
- If you get alternative site blood sugar results that do not agree with how you feel
- During illness or times of stress.
- If you will be driving a car or operating machinerv.

WARNING: Potential Biohazard

 The MICROLET NEXT lancing device is intended for single patient use only. It should not be used for assisted blood draws by health care

providers or at health care provision sites, and should never be shared with anyone else, even a family member, due to risk of infection.

- The lancing device is for use only on a single patient. Disinfect reusable components according to these instructions between each use
- Do not use on more than one patient. Improper use of blood lancets or the lancing device can increase the risk of inadvertent transmission of bloodborne pathogens, particularly in settings where multiple patients are tested. The cleaning and disinfection instructions for this device are intended only to reduce the risk of local use site infection: they cannot render this device safe for use for more than one patient.
- The lancet is for use only on a single patient. Discard the entire lancet after use. Used lancets are not sterile. Use a new sterile lancet each time vou test.
- Used lancets must be safely discarded after a single use. Always dispose of used test strips and lancets as medical waste or as advised by vour health care professional.
- Always wash your hands with soap and water and dry them well before and after testing or handling the meter, lancing device, or test strips.



Prepare the Lancing Device INTENDED USE: he MICROLET®NEXT lancing device is WARNING intended for collecting Keep out of reach of children. This kit contains capillary blood from small parts that could cause suffocation if the fingertip or palm for accidentally swallowed. blood glucose testing Lancet remains sterile until the protective cap is or other testing utilizing removed. Do not use if the lancet is damaged small amounts of blood or opened unintentionally. Dispose as medical The MICROLET NEXT lancing device is for waste to prevent any possible injury. If an individual is under the age of 18, the single-patient use. The MICROLET NEXT lancing lancing device and lancet must be used by an

- adult or under the supervision of the adult. . Wash hands and the puncture site with soap and warm water. Rinse and dry well. 2. Unlock the endcap. Twist the
- white band so the white raised mark and black indented mark line up as shown




"Click" to unlock

Pull off the endcap.

Loosen the round protective cap on the lancet by rotating the cap 1/4 turn, but **do not remove** the cap.



5. Insert the lancet firmly into the lancing device. **Push** until the lancet comes to a complete stop.



6. Carefully twist off the round protective cap. Save the cap for disposing of the used lancet.



for collecting fingertip capillary blood from persons age two years and older. The alternative site testing endcap is intended for collecting capillary blood from the palm in adults only

May not be included in all kits. Available on request from Customer Service


Get the Blood Drop - Fingertip

FINGERTIP TESTING

- Use the **black/white** adjustable endcap. Do not use the clear AST endcap.
- Attach the **black/white** adjustable endcap. Lock the endcap. Twist the white band until you
- hear or feel a click and the white and black raise marks line up as shown.



 Rotate the endcap dial to adjust the puncture depth (1.25 mm to 2.75 mm). The difference between each setting is 0.4 mm or less. The amount of pressure applied to the puncture site also affects puncture depth.



Setting	Puncture Depth (mm)	Setting	Puncture Depth (mm)
1	1.25	4	2.38
2	1.65	5	2.75
3	2.00		

Before using the lancing device for the first time on children, ask vour healthcare professional. Start with the shallowest puncture depth. If the puncture depth is too deep, the puncture can injure children. 4. Pull back the re-set handle to prepare the lancing device. The release button is raised when the lancing device is primed.

Press the lancing device firmly against the side of the finger and press the release button UNTIL THE DEVICE CLICKS.



6. Stroke your hand and finger toward the puncture site to form a drop of blood. Do not squeeze around the puncture site.



7. Test immediately after you have formed a small round blood drop, as shown.



Wipe the puncture site with a clean, dry tissue and **maintain pressure** on the site until the bleeding stops.

Eject / Dispose of the Lancet

WARNING: Potential Biohazard Used lancets must be safely discarded after a single use. Always dispose of used lancets as medical waste or as advised by your health care professional.

Your MICROLET NEXT lancing device has a lancet ejector. Do not use your fingers to remove the lancet from the lancing device.





. Carefully remove the endcap from the lancing device. Twist the white band of the endcap so the white raised mark and the black indented mark line up. Pull off the endcap.



"Click" to unlock

2. Place the round protective cap that you saved onto a flat surface. With the used lancet still in the lancing device, push the lancet into the center of the cap.





Continued on reverse side